



Preparation of cell extracts

 Toshiaki Tanaka*  Koji Moriya

Updated date: Jul 15, 2022

*For correspondence: ttanaka@bio.titech.ac.jp

 An abbreviated version of this protocol was published in Life Science Alliance in Feb 2022

Visualized procollagen Iα1 demonstrates the intracellular processing of propeptides

DOI: [10.26508/lsa.202101060](https://doi.org/10.26508/lsa.202101060)

Detailed protocol

Forty-eight hours after transfection, cells in 35 mm dish were washed twice with PBS and lysed with 50 µL of EBC lysis buffer (50 mM Tris-HCl, pH 8.0, 150 mM NaCl, 1 mM EDTA, and 0.5% NP-40) supplemented with an inhibitor mix (1 ng/ml of aprotinin, 100 mM β-glycerophosphate, 1 mM NaF, 1 mM Na3VO4, 10 µg/ml of leupeptin, 10 µg/ml of pepstatin A, and 1 mM phenylmethylsulfonyl fluoride) or 50 µL of NativePAGE sample buffer (Thermo Fisher Scientific) supplemented with 10% n-dodecyl-β-D-maltoside and the inhibitor mix described above. The cell extracts were cleared by centrifugation at 15,000 rpm for 15 min at 4 °C. The protein concentration of the precleared extracts was measured with Pierce BCA protein assay kit (Thermo Fisher Scientific).

How to cite: (Readers should cite both the Bio-protocol preprint and the original research article where this protocol was used)

1. Tanaka, T. and Moriya, K. (2022). Preparation of cell extracts. Bio-protocol Preprint. bio-protocol.org/prep1797.
2. Tanaka, T., Moriya, K., Tsunenaga, M., Yanagawa, T., Morita, H., Minowa, T., Tagawa, Y., Hanagata, N., Inagaki, Y. and Ikoma, T. (2022). Visualized procollagen Iα1 demonstrates the intracellular processing of propeptides. Life Science Alliance 5(5). DOI: [10.26508/lsa.202101060](https://doi.org/10.26508/lsa.202101060)

Copyright: Content may be subjected to copyright.